

# **Generation of $^{32}\text{P}$ -labeled MEK, ERK and VEGF-R protein standards for 2D gels**

*Jon Johansen, Matt Hoelter & Nancy Kendrick\**

Kendrick Labs Inc, Madison, WI


[www.kendricklabs.com](http://www.kendricklabs.com)

# Outline

- ❖ Protein kinase genes in the human genome
- ❖  $^{32}\text{P}$  labeling of ERK with MEK, a recombinant serine kinase
- ❖  $^{32}\text{P}$  labeling of VEGF-R, a recombinant tyrosine kinase
- ❖ Conclusions

## About Protein Kinases:

- ❖ The human genome contains about 500 kinase genes:
  - ~410 are serine/threonine kinases
  - ~90 are tyrosine kinases
  - most are not characterized
- ❖ Rough estimates suggest that 30-50% of proteins in mammalian cell lysates are serine/threonine phosphorylated. Tyrosine phosphorylation is rare and associated with mitotic events.

- 
- ❖ Protein standards known to be phosphorylated at specific residues would be useful. But none are commercially available. (Mass spec. labs use phosphorylated *peptides*)
  - ❖ So we decided to make our own phosphorylated standards using recombinant kinases and substrates.

## New company called ProQinase!

❖ specializes in recombinant protein kinases

❖ based in Freiburg, Germany but has US distributor



The screenshot shows a Windows Internet Explorer browser window displaying the ProQinase website. The address bar shows the URL <http://www.proqinase.com/>. The website features a logo on the left consisting of a stylized grey 'P' with a red dot. To the right of the logo, the text 'ProQinase' is written in a large, bold, red font, with 'TOOLS & TESTS' in a smaller, grey font below it. A navigation menu below the logo contains five items: 'Home', 'Company', 'Products & Services', 'Science', and 'Contact', each accompanied by a small image. The main content area has a heading 'Focus on Protein Kinases for Drug Development' followed by a paragraph describing ProQinase GmbH as a contract research organization (CRO) providing a 'Protein Kinase Technology Platform' for preclinical drug development. A second paragraph details the company's services, including subcutaneous xenograft and orthotopic tumor models, and a clinical biomarker assay service. On the right side of the page, there is a large image of a man in a lab coat looking at a microplate. Overlaid on this image are three promotional text blocks: 'New! Sets of protein kinases available Kinase Sampler and Kinase Sampler Plus', 'New Sophisticated Soft Agar Screening Services! Soft Agar Assay', and 'New Sophisticated Cellular Assay Services! Cell Migration Assay'.

protein kinase - kinase assay - angiogenesis assay - Windows Internet Explorer

<http://www.proqinase.com/>

ProQinase  
TOOLS & TESTS

Home Company Products & Services Science Contact

### Focus on Protein Kinases for Drug Development

ProQinase GmbH is a contract research organization (CRO) which provides a *Protein Kinase Technology Platform* for preclinical drug development of protein kinase inhibitors in oncology and other therapeutic areas. Currently 185 in-house produced recombinant protein kinases are offered for sale and 320 kinases are available for *in vitro* kinase assay services (HTS, selectivity profiling and IC<sub>50</sub> determination). Hits derived from *in vitro* kinase assays can be characterized at ProQinase in cellular phosphorylation assays and in a cellular angiogenesis assay.


For further evaluation of lead compounds ProQinase offers various *in vivo* assays including subcutaneous xenograft and orthotopic tumor models, and an unique angiogenesis assay. ProQinase's drug development service portfolio is completed by a clinical biomarker assay service for determination of different biomarker proteins in serum or plasma samples derived from clinical trials.

**New!**  
Sets of protein kinases available  
Kinase Sampler and Kinase Sampler Plus

**New Sophisticated Soft Agar Screening Services!**  
Soft Agar Assay

**New Sophisticated Cellular Assay Services!**  
Cell Migration Assay

## Recombinant MEK kinase and ERK substrate were purchased from ProQinase.



MEK is a phosphoserine kinase involved in signal transduction.

## The phosphorylation reaction was carried out in vitro


MEK-1 kinase (MW 43,600 from Sf9 insect cells)

ERK-2 substrate (MW 45,500 from *E. coli*)

ProKinase Standard Assay Buffer

Either cold ATP or  $^{32}\text{P}$ -labeled ATP

The reaction was allowed to proceed 40 min at 30° C, stopped with  $\text{H}_3\text{P}_4$ , diluted with urea buffer and aliquots run on 2D gels.




The kinase mix was resolved by 2D Electrophoresis and Western blotting with anti-PTyr or anti-PSer/-PThr antibodies.

- ❖ A 2D gel is run, and the proteins transferred to PVDF. The membrane is stained with Coomassie blue and scanned, shaken overnight with an antibody, treated with ECL or ECL Plus and exposed to x-ray film
- ❖ The antibody is either:
  - A generic monoclonal antibody, PY20, against phosphotyrosine or
  - Qiagen monoclonal antibodies, Q5 against phosphoserine, Q7 against phosphothreonine

More method details at:  
[www.kendricklabs.com](http://www.kendricklabs.com)









## Silver stain versus Western blot



The main ERK isoform moves in the acidic direction during IEF, suggesting phosphorylation is occurring when the MEK kinase is present, as expected. Surprisingly, the anti-phosphoserine ab lights up the MEK kinase as well as both forms of the Erk substrate.


## Coomassie and WB versus $^{32}\text{P}$ label

MEK + ERK + $^{32}\text{PO}_4$	CB-PVDF	
ERK + $^{32}\text{PO}_4$	CB-PVDF	
MEK + ERK + $^{32}\text{PO}_4$	x-ray film	 43
ERK + $^{32}\text{PO}_4$	x-ray film	 43
MEK + ERK + $^{32}\text{PO}_4$	ECL-WB	 43
ERK + $^{32}\text{PO}_4$	ECL-WB	 43

Two 2D gels were run: MEK kinase + ERK substrate +  $^{32}\text{P}$ , and substrate alone +  $^{32}\text{P}$ . The blots were exposed to x-ray film and then P-Ser WB was carried out.

*Conclusions: The MEK kinase is clearly phosphorylated. The unphosphorylated ERK substrate lights up with the anti-phosphoserine antibody - non-specific binding, unless ERK is phosphorylated by the insect cell line.*


## Carbamylated proteins provide evidence that 2D gels show single charge changes



NL Anderson and BJ Hickman, *Anal Biochem*, 1979; 93: 312.

# Quantitative Analysis

	Spot #	CB Vol	Film Vol	Film/CB Specific Activity	Ratio to most basic
MEK	1	0.5	3	7	4
	2	1.1	32	28	14
	3	6	110	17	9
	4	17	156	9	5
	5	21	100	5	3
	6	10	18	2	1
ERK	7	1.1	14	12	1
	8	10	232	23	2
	9	3	78	28	3
	10	1.1	12	11	1
	11	0.4	3	10	1



$^{32}\text{P}$ -labeled MEK/ERK was run on a Coomassie-stained 2D gel and exposed to x-ray film. Both gel and film were quantitatively scanned and analyzed with Progenesis SameSpots software from Nonlinear Dynamics. But results didn't make sense in terms of phosphates/isoform.


First assumption, that each aa site reaches 100% phosphorylation before the next one is filled, is naive.

MEK	% AA phosphorylation			
	Condition 1	0%	0%	100%
10 molecules	aa 12	aa 48	aa 96	aa 192
1			P-Ser	
2			P-Ser	
3			P-Ser	
4			P-Ser	
5			P-Ser	
6			P-Ser	
7			P-Ser	
8			P-Ser	
9			P-Ser	
10			P-Ser	

In reality... variable phosphorylation of each aa site. We'll have to send the individual spots to our collaborator at Columbia University Protein Core for analysis to determine phosphates per isoform.

MEK	% AA phosphorylation			
	Condition 4	10%	30%	60%
10 molecules	aa 12	aa 48	aa 96	aa 192
1	P-Ser	P-Ser	P-Ser	P-Ser
2	20%	P-Ser	P-Ser	P-Ser
3		P-Ser	P-Ser	P-Ser
4		30%	P-Ser	P-Ser
5			P-Ser	P-Ser
6			P-Ser	P-Ser
7			30%	P-Ser
8				P-Ser
9				P-Ser
10				10%

# Vascular Endothelial Growth Factor Receptor (VEGF)




VEGF is a transmembrane receptor

- ❖ VEGF stimulates the formation of a new blood supply for developing tumors via its receptor. Inhibitors of VEGF-R (Avastin/Genentech) have been approved to treat a variety of cancers.
- ❖ We bought VEGF-R from ProQinase,  $^{32}\text{P}$  labeled assuming self-phosphorylation, TCA ppt'd to concentrate, ran a 2D gel, transferred to PVDF, exposed 4 days and then did an anti-PTyr Western blot.

## The PY20 antibody is *very* sensitive!



VEGF-R <sup>32</sup>P Autorad: 4 day x-ray film exposure to show <sup>32</sup>P labeled protein. Nothing was visible by Coomassie blue on the PVDF.



VEGF-R Western blot: 30 sec film exposure after anti-PTyr antibody exposure and ECL treatment.

We'll buy and label a non-receptor tyrosine kinase to use as an internal standard.

# Conclusions

- ❖ Generating phosphorylated protein standards using recombinant proteins from ProKinase is straightforward.
- ❖ Phosphorylation shifts protein isoforms on 2D gels in the acidic direction as expected. However, the extent of the shift doesn't seem proportional to the amount of  $^{32}\text{P}$  incorporated. Mass spectrometry will be necessary to determine # of phosphates/isoform.
- ❖ The final standards should have many uses.